

R J Hill Laboratories Limited 1 Clyde Street Private Bag 3205 Hamilton 3240, New Zealand

Tel +64 7 858 2000 Fax +64 7 858 2001 Email mail@hill-labs.co.nz Web www.hill-labs.co.nz

ANALYSIS REPORT

Page 1 of 2

SPv1

Client: Hawkes Bay Regional Council

Contact: V Lyon

C/- Hawkes Bay Regional Council

Private Bag 6006 NAPIER 4142 Lab No:1422201Date Registered:06-May-2015Date Reported:11-May-2015

Quote No: 66907 Order No: N37006 Client Reference: Whangawehi Add. Client Ref: 312-302 V Lyon

Sample Type: Aqueous						
Sa	ample Name:	54372 - Whangawehi Strm at Pat O'Brians-3304	54373- Mangatupae Strm at Pat O'Brians-3303	54374 - Whangawehi at George Ormonds -3301	54375 - Coops Crk US W hangawhi - 3306	54376 - Reserve Stream - Trib of Whangawehi - 3307
	Lab Number:	1422201.1	1422201.2	1422201.3	1422201.4	1422201.5
Individual Tests						
Escherichia coli	cfu / 100mL	5 ^{#1}	50 #1	3 #1	40 #1	90 #1
HBRC Standard River Trace Nu	trients					
Total Suspended Solids	g/m³	0.7	< 0.5	0.6	0.7	1.6
Total Nitrogen	g/m³	0.21	0.18	0.23	0.13	0.20
Total Ammoniacal-N	g/m³	< 0.005	< 0.005	< 0.005	< 0.005	0.021
Nitrite-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.0023
Nitrate-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.056
Nitrate-N + Nitrite-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.058
Total Kjeldahl Nitrogen (TKN)	g/m³	0.21	0.18	0.23	0.13	0.14
Dissolved Reactive Phosphorus	g/m³	0.030	0.0075	0.034	0.020	0.073
Total Phosphorus	g/m³	0.038	0.017	0.042	0.025	0.107
Sa	ample Name:	54377 - Whangawehi US Reserve Confl - 3308	54378 - Whangawehi DS Cattleyards - 3309			
	Lab Number:	1422201.6	1422201.7			
Individual Tests						
Escherichia coli	cfu / 100mL	120 #1	60 #1	-	-	-
HBRC Standard River Trace Nu	trients					
Total Suspended Solids	g/m³	0.8	1.3	-	-	-
Total Nitrogen	g/m³	0.43	0.43	-	-	-
Total Ammoniacal-N	g/m³	0.010	< 0.005	-	-	-
Nitrite-N	g/m³	0.0026	< 0.0010	-	-	-
Nitrate-N	g/m³	0.24	0.0097	-	-	-
Nitrate-N + Nitrite-N	g/m³	0.25	0.0103	-	-	-
Total Kjeldahl Nitrogen (TKN)	g/m³	0.19	0.42	-	-	-
Dissolved Reactive Phosphorus	g/m³	0.071	0.024	-	-	-
Total Phosphorus	g/m³	0.091	0.045	-	-	-

Analyst's Comments

Please interpret these microbiological results with caution as the age of the sample is unknown. Samples must be less than 24 hours old at the time of testing. Please record both the time and date of sampling on your submission form and the sample bottles.

^{#1} Statistically estimated count based on the theoretical countable range for the stated method.





SUMMARY OF METHODS

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively clean matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis.

Sample Type: Aqueous	Mathad Decembrian	Default Detection 1 limit	Commis No
Test	Method Description		Sample No
HBRC Standard River Trace Nutrients		0.0010 - 0.5 g/m ³	1-7
Filtration, Unpreserved	Sample filtration through 0.45µm membrane filter.	-	1-7
Total Kjeldahl Digestion	Sulphuric acid digestion with copper sulphate catalyst.	-	1-7
Total Phosphorus Digestion	Acid persulphate digestion.	-	1-7
Total Suspended Solids	Filtration of a 2L sample using Whatman 934 AH, Advantec GC-50 or equivalent filters (nominal pore size 1.2 - 1.5µm), gravimetric determination. APHA 2540 D 22 nd ed. 2012.	0.5 g/m ³	1-7
Total Nitrogen	Calculation: TKN + Nitrate-N + Nitrite-N. Please note: The Default Detection Limit of 0.05 g/m³ is only attainable when the TKN has been determined using a trace method utilising duplicate analyses. In cases where the Detection Limit for TKN is 0.10 g/m³, the Default Detection Limit for Total Nitrogen will be 0.11 g/m³.	0.05 g/m³	1-7
Total Ammoniacal-N Trace	Phenol/hypochlorite colorimetry. Flow injection analyser. (NH4-N = NH4+-N + NH3-N). APHA 4500-NH ₃ H 22 nd ed. 2012.	0.005 g/m ³	1-7
Nitrite-N Trace	Automated Azo dye colorimetry, Flow injection analyser. APHA 4500-NO ₃ ·I 22 nd ed. 2012.	0.0010 g/m ³	1-7
Nitrate-N	Calculation: (Nitrate-N + Nitrite-N) - NO2N. In-House.	0.0010 g/m ³	1-7
Nitrate-N + Nitrite-N Trace	Total oxidised nitrogen. Automated cadmium reduction, flow injection analyser. APHA 4500-NO ₃ ·I 22 nd ed. 2012.	0.0010 g/m ³	1-7
Total Kjeldahl Nitrogen (TKN)	Total Kjeldahl digestion, phenol/hypochlorite colorimetry. Discrete Analyser. APHA 4500-N _{org} D. (modified) 4500 NH ₃ F (modified) 22 nd ed. 2012.	0.10 g/m ³	1-7
Dissolved Reactive Phosphorus (trace)	Filtered sample. Molybdenum blue colorimetry. Flow injection analyser. APHA 4500-P G 22 nd ed. 2012.	0.0010 g/m ³	1-7
Total Phosphorus	Total phosphorus digestion, ascorbic acid colorimetry. Discrete Analyser. APHA 4500-P B & E (modified from manual analysis) 22^{nd} ed. 2012. Also modified to include the use of a reductant to eliminate interference from arsenic present in the sample. NWASCA, Water & soil Miscellaneous Publication No. 38, 1982.	0.004 g/m ³	1-7
Escherichia coli	Membrane filtration, Count on mFC agar, Incubated at 44.5°C for 22 hours, MUG Confirmation. Analysed at Hill Laboratories - Microbiology; 1 Clow Place, Hamilton. APHA 9222 G, 22 nd ed. 2012.	1 cfu / 100mL	1-7

These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Samples are held at the laboratory after reporting for a length of time depending on the preservation used and the stability of the analytes being tested. Once the storage period is completed the samples are discarded unless otherwise advised by the client.

This report must not be reproduced, except in full, without the written consent of the signatory.

Ara Heron BSc (Tech)

Client Services Manager - Environmental Division