

R J Hill Laboratories Limited 1 Clyde Street Private Bag 3205 Hamilton 3240, New Zealand

Tel +64 7 858 2000 Fax +64 7 858 2001 Email mail@hill-labs.co.nz Web www.hill-labs.co.nz

ANALYSIS REPORT

Page 1 of 2

SPv1

Client: H

Hawkes Bay Regional Council

Contact: V Lyon

C/- Hawkes Bay Regional Council

Private Bag 6006 NAPIER 4142 Lab No: 1519254

Date Registered: 23-Dec-2015

Date Reported: 12-Jan-2016

Guote No: 69783

Order No: N43149

Client Reference: Whangawehi

Add. Client Ref: 312-302 Submitted By: V Lyon

Sample Type: Aqueous						
Sar	mple Name:	58506 - W hangawehi Strm at Pat O'Brians-3304	58507 - Mangatupae Strm at Pat O'Brians-3303	58508 - Whangawehi at George Ormonds -3301	58509 - Coops - Trib of Whangawhi - 3306	58510 - Reserve Stream - Trib of Whangawehi - 3307
La	ab Number:	1519254.1	1519254.2	1519254.3	1519254.4	1519254.5
Individual Tests				,		
Escherichia coli	cfu / 100mL	630 #2	74 #2	80 #2	160 #2	60 #2
HBRC Standard River						
Volatile Suspended Solids	g/m ³	0.6 #1	0.6	0.8	< 0.5	< 0.5
Total Suspended Solids	g/m³	0.5 #1	1.4	1.0	0.8	1.0
Total Nitrogen	g/m³	0.27	0.24	0.30	0.21	0.36
Total Ammoniacal-N	g/m³	< 0.005	< 0.005	< 0.005	< 0.005	< 0.005
Nitrite-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.0012
Nitrate-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.048
Nitrate-N + Nitrite-N	g/m³	< 0.0010	< 0.0010	< 0.0010	< 0.0010	0.049
Total Kjeldahl Nitrogen (TKN)	g/m³	0.27	0.24	0.30	0.21	0.31
Dissolved Reactive Phosphorus	g/m³	0.036	0.0192	0.020	0.033	0.056
Total Phosphorus	g/m³	0.046	0.030	0.038	0.042	0.064
Sar	nple Name:	58511 -	58512 -			
		Whangawehi US Reserve Confl - 3308	Whangawehi DS Cattleyards - 3309			
1:	ab Number:	1519254.6	1519254.7			
Individual Tests				I.		I.
Escherichia coli	cfu / 100mL	40	350	-	-	-
HBRC Standard River						
Volatile Suspended Solids	g/m³	< 0.5	< 0.5	-	-	-
Total Suspended Solids	g/m³	0.9	0.6	-	-	-
Total Nitrogen	g/m ³	0.17	0.40	-	-	-
Total Ammoniacal-N	g/m ³	< 0.005	< 0.005	-	-	-
Nitrite-N	g/m ³	0.0015	< 0.0010	-	-	-
Nitrate-N	g/m ³	0.024	< 0.0010	-	-	-
Nitrate-N + Nitrite-N	g/m ³	0.025	< 0.0010	-	-	-
Total Kjeldahl Nitrogen (TKN)	g/m ³	0.15	0.40	-	-	-
Dissolved Reactive Phosphorus	g/m ³	0.094	0.036	-	-	-
Total Phosphorus	g/m ³	0.105	0.067	-	-	-

Analyst's Comments

- ^{#1} It has been noted that the result for Volatile Suspended Solids was greater than that for Total Suspended Solids, but within the analytical variation of these methods.
- #2 Statistically estimated count based on the theoretical countable range for the stated method.



SUMMARY OF METHODS

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively clean matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis.

Sample Type: Aqueous						
Test	Method Description	Default Detection Limit	Sample No			
HBRC Standard River		0.0010 - 0.5 g/m ³	1-7			
Filtration, Unpreserved	Sample filtration through 0.45µm membrane filter.	-	1-7			
Total Kjeldahl Digestion	Sulphuric acid digestion with copper sulphate catalyst.	-	1-7			
Total Phosphorus Digestion	Acid persulphate digestion.	-	1-7			
Volatile Suspended Solids	Filtration (GF/C, 1.2 μm). Ashing 550°C, 30 min. Gravimetric. APHA 2540 E 22 nd ed. 2012.	0.5 g/m ³	1-7			
Total Suspended Solids	Filtration of a 2L sample using Whatman 934 AH, Advantec GC-50 or equivalent filters (nominal pore size 1.2 - 1.5µm), gravimetric determination. APHA 2540 D 22 nd ed. 2012.	0.5 g/m ³	1-7			
Total Nitrogen	Calculation: TKN + Nitrate-N + Nitrite-N. Please note: The Default Detection Limit of 0.05 g/m³ is only attainable when the TKN has been determined using a trace method utilising duplicate analyses. In cases where the Detection Limit for TKN is 0.10 g/m³, the Default Detection Limit for Total Nitrogen will be 0.11 g/m³.	0.05 g/m³	1-7			
Total Ammoniacal-N Trace	Phenol/hypochlorite colorimetry. Flow injection analyser. (NH4-N = NH4+-N + NH3-N). APHA 4500-NH3 H 22 nd ed. 2012.	0.005 g/m ³	1-7			
Nitrite-N Trace	Automated Azo dye colorimetry, Flow injection analyser. APHA 4500-NO ₃ -I 22 nd ed. 2012 (modified).	0.0010 g/m ³	1-7			
Nitrate-N	Calculation: (Nitrate-N + Nitrite-N) - NO2N. In-House.	0.0010 g/m ³	1-7			
Nitrate-N + Nitrite-N Trace	Total oxidised nitrogen. Automated cadmium reduction, flow injection analyser. APHA 4500-NO ₃ · I 22 nd ed. 2012 (modified).	0.0010 g/m³	1-7			
Total Kjeldahl Nitrogen (TKN)	Total Kjeldahl digestion, phenol/hypochlorite colorimetry. Discrete Analyser. APHA 4500-N _{org} D. (modified) 4500 NH ₃ F (modified) 22 nd ed. 2012.	0.10 g/m ³	1-7			
Dissolved Reactive Phosphorus (trace)	Filtered sample. Molybdenum blue colorimetry. Flow injection analyser. APHA 4500-P G 22 nd ed. 2012.	0.0010 g/m³	1-7			
Total Phosphorus	Total phosphorus digestion, ascorbic acid colorimetry. Discrete Analyser. APHA 4500-P B & E (modified from manual analysis) 22 nd ed. 2012. Also modified to include the use of a reductant to eliminate interference from arsenic present in the sample. NWASCA, Water & soil Miscellaneous Publication No. 38, 1982.	0.004 g/m ³	1-7			
Escherichia coli	Membrane filtration, Count on mFC agar, Incubated at 44.5°C for 22 hours, MUG Confirmation. Analysed at Hill Laboratories - Microbiology; 1 Clow Place, Hamilton. APHA 9222 G, 22 nd ed. 2012.	1 cfu / 100mL	1-7			

These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Samples are held at the laboratory after reporting for a length of time depending on the preservation used and the stability of the analytes being tested. Once the storage period is completed the samples are discarded unless otherwise advised by the client.

This report must not be reproduced, except in full, without the written consent of the signatory.

Ara Heron BSc (Tech)

Client Services Manager - Environmental Division