

Private Bag 3205 Hamilton 3240 New Zealand

E mail@hill-labs.co.nz W www.hill-laboratories.com

NALYSIS REPORT

Page 1 of 2

SPv1

Client: Hawkes Bay Regional Council

Contact: V Lyon

C/- Hawkes Bay Regional Council

Private Bag 6006 Napier 4142

1838099 Lab No: **Date Received:** 06-Sep-2017 **Date Reported:** 20-Sep-2017 **Quote No:** 78490 **Order No:** N46928 Whangawehi Client Reference: 312-302 Add. Client Ref: **Submitted By:** V Lyon

Sample Type: Aqueous								
S	Sample Name:	63499 - Coops -	63500 - Reserve	63501 -				
		Trib of	Stream - Trib of	Whangawehi US				
		Whangawhi -	Whangawehi -	Reserve Confl -				
		3306	3307	3308				
		05-Sep-2017	05-Sep-2017	05-Sep-2017 1:00				
		12:30 pm	12:45 pm	pm				
	Lab Number:	1838099.1	1838099.2	1838099.3				
Faecal Coliforms and E. coli pro	ofile							
Faecal Coliforms	cfu / 100mL	650	450	390	-	-		
Escherichia coli	cfu / 100mL	630	330	370	-	-		
HBRC Standard River								
Volatile Suspended Solids	g/m³	2.7	2.2	1.2	-	-		
Total Suspended Solids	g/m³	32	32	9.6	-	-		
Total Nitrogen	g/m³	1.03	0.57	1.28	-	-		
Total Ammoniacal-N	g/m³	0.029	0.019	< 0.005	-	-		
Nitrite-N	g/m³	0.0051	0.0011	0.0034	-	-		
Nitrate-N	g/m³	0.69	0.24	0.85	-	-		
Nitrate-N + Nitrite-N	g/m³	0.69	0.24	0.85	-	-		
Total Kjeldahl Nitrogen (TKN)	g/m³	0.34	0.33	0.43	-	-		
Dissolved Reactive Phosphorus	g/m³	0.072	0.041	0.048	-	-		
Total Phosphorus	g/m³	0.083	0.054	0.070	-	-		

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively clean matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis

Sample Type: Aqueous						
Test	Method Description	Default Detection Limit	Sample No			
Individual Tests						
Filtration, Unpreserved	Sample filtration through 0.45µm membrane filter. Analysis performed at 1 Clyde Street, Hamilton.		1-3			
Total Kjeldahl Digestion	Sulphuric acid digestion with copper sulphate catalyst.	-	1-3			
Total Phosphorus Digestion	Acid persulphate digestion.	-	1-3			
Volatile Suspended Solids	Filtration (GF/C, 1.2 μ m). Ashing 550°C, 30 min. Gravimetric. APHA 2540 E 22 nd ed. 2012.	0.5 g/m ³	1-3			
Total Suspended Solids	Filtration of a 2L sample using Whatman 934 AH, Advantec GC-50 or equivalent filters (nominal pore size 1.2 - 1.5µm), gravimetric determination. APHA 2540 D 22 nd ed. 2012.	0.5 g/m ³	1-3			
Total Nitrogen	Calculation: TKN + Nitrate-N + Nitrite-N. Please note: The Default Detection Limit of 0.05 g/m³ is only attainable when the TKN has been determined using a trace method utilising duplicate analyses. In cases where the Detection Limit for TKN is 0.10 g/m³, the Default Detection Limit for Total Nitrogen will be 0.11 g/m³.	0.05 g/m³	1-3			



This Laboratory is accredited by International Accreditation New Zealand (IANZ), which represents New Zealand in the International Laboratory Accreditation Cooperation (ILAC). Through the ILAC Mutual Recognition Arrangement (ILAC-MRA) this accreditation is internationally recognised.

Sample Type: Aqueous							
Test	Method Description	Default Detection Limit	Sample No				
Total Ammoniacal-N Trace	Phenol/hypochlorite colorimetry. Flow injection analyser. (NH4-N + NH3-N). Analysis performed at 1 Clyde Street, Hamilton. APHA 4500-NH ₃ H 22 nd ed. 2012.		3				
Total Ammoniacal-N Trace	Phenol/hypochlorite colorimetry. Flow injection analyser. (NH4-N = NH4+-N + NH3-N). APHA 4500-NH3 H 22 nd ed. 2012.	0.005 g/m ³	1-2				
Nitrite-N Trace	Automated Azo dye colorimetry, Flow injection analyser. Analysis performed at 1 Clyde Street, Hamilton. APHA 4500-NO ₃ -I 22 nd ed. 2012 (modified).	0.0010 g/m ³	1-3				
Nitrate-N	Calculation: (Nitrate-N + Nitrite-N) - NO2N. In-House.	0.0010 g/m ³	1-3				
Nitrate-N + Nitrite-N Trace	Total oxidised nitrogen. Automated cadmium reduction, flow injection analyser. Analysis performed at 1 Clyde Street, Hamilton. APHA 4500-NO ₃ · I 22 nd ed. 2012 (modified).	0.0010 g/m ³	1-3				
Total Kjeldahl Nitrogen (TKN)	Total Kjeldahl digestion, phenol/hypochlorite colorimetry. Discrete Analyser. Analysis performed at 1 Clyde Street, Hamilton. APHA 4500-Norg D. (modified) 4500 NH ₃ F (modified) 22 nd ed. 2012.	0.10 g/m³	1-3				
Dissolved Reactive Phosphorus (trace)	Filtered sample. Molybdenum blue colorimetry. Flow injection analyser. Analysis performed at 1 Clyde Street, Hamilton. APHA 4500-P G 22 nd ed. 2012.	0.0010 g/m ³	1-3				
Total Phosphorus	Total phosphorus digestion, ascorbic acid colorimetry. Discrete Analyser. APHA 4500-P B & E (modified from manual analysis) 22 nd ed. 2012. Also modified to include the use of a reductant to eliminate interference from arsenic present in the sample. NWASCA, Water & soil Miscellaneous Publication No. 38, 1982.	0.004 g/m ³	1-3				
HBRC Standard River		0.0010 - 0.5 g/m ³	1-3				
Faecal Coliforms and E. coli profile		1					
Faecal Coliforms	Membrane Filtration, Count on mFC agar, Incubated at 44.5°C for 22 hours, Confirmation. Analysed at Hill Laboratories - Microbiology; 1 Clow Place, Hamilton. APHA 9222 D, 22 nd ed. 2012.	1 cfu / 100mL	1-3				
Escherichia coli	Membrane filtration, Count on mFC agar, Incubated at 44.5°C for 22 hours, MUG Confirmation. Analysed at Hill Laboratories - Microbiology; 1 Clow Place, Hamilton. APHA 9222 G, 22 nd ed. 2012.	1 cfu / 100mL	1-3				

These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Samples are held at the laboratory after reporting for a length of time depending on the preservation used and the stability of the analytes being tested. Once the storage period is completed the samples are discarded unless otherwise advised by the client.

This report must not be reproduced, except in full, without the written consent of the signatory.

Ara Heron BSc (Tech)

Client Services Manager - Environmental